

PRODUCT DESCRIPTION

Cellaria’s WIT-P™ Culture Medium supports the expansion of human primary mammary epithelial cells for more than 20 passages without the growth arrest that is seen with other primary epithelial cell culture media. In addition, primary mammary epithelial cells grown in WIT-P Culture Medium exhibit a more normal luminal pattern of gene expression than cells grown in other epithelial culture media. Additional applications for WIT in tumor and normal cell biology can be found in peer-reviewed publications.



ORDERING INFORMATION

<i>PRODUCT NAME</i>	<i>CATALOG NUMBER</i>	<i>COMPONENTS</i>	<i>PART NUMBER</i>	<i>SIZE</i>
WIT-P	CM-0101	WIT Basal Medium	CM-0100-A1	500 mL
		WIT-P Supplement	CM-0101-P1	15 mL

PRODUCT HANDLING/DIRECTIONS FOR USE

Reconstitution:

To make WIT-P Complete Medium, please combine the following:

1. 500 mL WIT Basal Medium
2. 15 mL WIT-P Supplement (thawed in a 37°C water bath)

Storage & Stability:

Upon arrival, store basal medium at 4°C. Store supplement at -20°C. Medium and supplement are stable for a minimum of 3 months from date of receipt when stored as directed.